

medinov

99, rue de Gerland

69007 Lyon

Tél: 04 37 28 08 14

www.medinov.fr

Fournisseur Agréé

N° IDENTIFICATION 69.2.65115.1

UPLIFT POWER SEAT

The electric lifting cushion that turns an armchair or sofa into a lift chair.

The Uplift Power Seat is an electric lifting cushion that gently seats and raises users from an armchair or sofa. It provides safe, controlled support and lifts 100% of the user's weight.

Power Seat is light weight (12 lb/5.4 kg), has a built-in carry handle, and is easily transported from room to room.

BENEFITS

- Use on existing furniture
- Lifts users weighing up to 300 lb (136 kg)
- Works on sofas, armchairs and recliners
- Costs considerably less than a lift chair
- Transports easily
- Provides a flat, stable lifting surface (will not pitch users forward)

WHO CAN BENEFIT

People who have lost sufficient upper or lower body strength to get up unassisted. Those who have lost mobility as a result of:

- Arthritis
- Parkinson's
- Muscular dystrophy
- ALS (Lou Gerhig's disease)
- Post operative conditions
- Degenerative joint disease



UPLIFT
POWER SEAT™
Bringing quality to life.

SPECIFICATIONS

Base Dimensions Depth: 17"/43 cm
Width: 14"/36 cm

Seat Dimensions Depth: 19"/48 cm
Width: 16"/41 cm

Seat Height Open: 10"/25 cm
Flat: 2"/5 cm

Weight Unit: 12 lb/5.4 kg.
User: Up to 300 lb/
136 kg

Power Supply 120 Vac 60Hz 0.6A
CSA & C/US Approved

Cord Length 10 ft/3 m



To rise, pull up on the Power Lever, and then place your hand(s) on the arm(s) of the chair or sofa while you are being lifted.

EASY CARE
Waterproof incontinent cover is easily removed and machine washed. Water resistant coating on foam surface can be wiped down with a damp cloth.



The Uplift Power Seat will stop automatically when it reaches its maximum height, or you can stop it by gently pushing down on the Power Lever. Flexible seat curves to flattened seating surface at every lift level, offering a stable and safe lift. Patent # 5,316,370 & 6,702,383. European and international patents pending.

SAFE

Power Lever
Flexible seat curves to flattened seating surface at every lift level, offering a stable and safe lift. Patent # 5,316,370 & 6,702,383. European and international patents pending.



Now you can leave the armchair or sofa. The Uplift Power Seat is ready to help seat you when you return to your comfortable sofa.

COMFORTABLE

Cushion is made of a high density, slow-release memory foam that helps protect against pressure sores and conforms to the shape of the body - providing superior comfort.

MEMORY FOAM

SECURE

Base has a non-slip material bottom that prevents the cushion from shifting.

EASY TO USE

Simply pull up on one of the two Power Levers to rise; push down to lower.



DURABLE

Durable construction accommodates users weighing up to 300 lb (136 kg) and is warranted for one year when used as recommended on armchairs or sofas (with arms).

POWER LEVERS

Two detachable Power Levers allow for use of either left or right hand.

TRANSPORTABLE

Built-in handle makes it easy to carry Uplift Power Seat from place to place.



UPLIFT
POWER SEAT
Bringing quality to life.

Manufactured for Uplift Technologies Inc.

medinov
99, rue de Gerland
69007 Lyon
Tél: 04 37 28 08 14
www.medinov.fr
Fournisseur Agréé
N° IDENTIFICATION 69.2.65115.1